# Aether MyoSense

## Instruction For Use

For Qualified Personnel

Version 2.0 / 2026.06.08

EN / PL / UA

# Aether Myo Sense

## Instruction For Use

### Table Of Contents

Aether MyoSense .....	2
Intended Use .....	2
Technical Specification .....	4
Scope of delivery .....	5
Features.....	5
Compatibility.....	6
Set up in Prosthesis.....	7
Configuration and Adjustments.....	8
Maintenance.....	9
Reliability and Safety .....	9
Warnings.....	10
Troubleshooting .....	11
Regulatory Compliance .....	13
Symbols .....	13
EMC Labelling Information.....	16
Reporting.....	17

We highly recommend reading this document carefully before using the device.

This document contains information about the proper and safe use of the device.

We highly recommend keeping this document for the whole period of using the device.

If you have any questions about the product, please contact us at [info@aetherbiomedical.com](mailto:info@aetherbiomedical.com).

## **Aether MyoSense**

Aether MyoSense is a dual-channel electromyography (EMG) signal acquisition and processing system designed to enable precise and intuitive control of upper-limb prosthetic devices. The system captures and processes muscle signals, translating voluntary muscle contractions into functional control signals.

Key features of the system include Bluetooth Low Energy (BLE) communication for configuration, real-time skin contact monitoring, and adjustable gain and control signal filtering options, ensuring that both clinicians and patients can optimize the device for individual needs.

Read this document carefully before fitting.

## **Intended Use**

Aether MyoSense is designed to amplify the signal received from the patient's skeletal muscles, utilizing surface electromyography. Aether MyoSense is an external device intended to be used with other appropriate upper limb components to form a complete arm prosthesis (class I device), to be fitted only by qualified and certified clinicians.

Enables the patient to control the end device (prosthesis) by delivering an electrical signal to the device proportional

to voluntary skeletal muscle contraction. Fitting should only be carried out by qualified and certified personnel.

## **Indications**

The patient must be able to understand usage and safety messages and put them into practice.

Amputation level below-elbow, above-elbow and shoulder and wrist disarticulations.

For unilateral or bilateral amputation.

Congenital upper limb deficiency.

## **Patient Population**

Aether MyoSense is recommended for all genders aged 14–75 years.

## **Contraindications**

The device is not suitable for patients who cannot generate sufficient voluntary muscle contractions for prosthetic control.

Not intended for use in industrial applications that could interfere with the surface electromyography readings.

The safety and effectiveness have not been established for pediatric or geriatric populations outside of the intended range.

The device is not suitable to be used in high electromagnetic interference environments exceeding norms considered as normal or home use.

The device cannot be used by people with cognitive deficits.

## **Storage and Handling**

Handle the device with care, avoiding excessive bending or pulling of wires and connectors. Avoid exposing it to sudden physical shocks or mechanical impacts (e.g.,

dropping). Avoid prolonged exposure to UV lightning, including sunlight.

## Safe usage

The final decision whether Aether MyoSense should be used belongs to qualified medical personnel.

Aether MyoSense is intended exclusively for use on one patient. Use of the product by another person is not approved by the manufacturer.


Fitting of the product to the patient's upper limb may be exclusively done by a qualified and certified clinician/prosthetist.

Aether MyoSense should not be used for the operation of motor vehicles, heavy equipment (e.g. construction machines), industrial machines, motor-driven equipment or any type of firearm or weapons.

Aether MyoSense is designed for mild to moderate activities.

## Technical Specification

Dimensions (main unit)	50 mm x 17.5 mm x 7 mm
Temperature range (use)	-5 °C to +45 °C
Temperature range (transport and storage)	-25 °C to +70 °C
Humidity range	5 % to 93 %
Atmospheric pressure range	700 hPa to 1060 hPa
Input voltage	6 to 13 VDC
Output voltage	0 to 5 VDC
Electrode leads length	25 cm
3-pin connectors length	15 cm

### Scope of delivery

Aether MyoSense [**AET-SEN-2NM**] is available in a single variation, dual-channel, with 15 cm 3-pin connector cables and 25 cm dome electrode leads.

### Features

#### Dual-channel electromyography (EMG) measurement

Aether MyoSense features two independent EMG measurement channels, translating the patient's skeletal muscle activity into two separate analogue control signals. Each channel consists of two differential electrode inputs. The fifth, reference (REF), electrode is used to ensure proper system performance.

## Connection to Aether Digital Platform (ADP)

Aether MyoSense is a Bluetooth (BLE) enabled device that can be connected to Aether Digital Platform (ADP) for adjustment and tuning. Establishing a BLE connection is possible up to 10 minutes after system power-up.

### Gain

The device supports digitally adjustable gain, allowing for accommodation of the control signal level to the patient's individual degree of voluntary contraction muscle activity.

### Smoothing

Aether MyoSense is fitted with an adjustable control signal filtering feature. This can be used to counteract sudden spikes and tremors in output signals, which prevents unintentional prosthetic device activation and eases the patient's ability to maintain a steady control signal, while keeping extremely low control latency.

### Skin contact quality

Aether MyoSense is capable of assessing the quality of contact between each channel's electrode pair and the patient's skin, expressed as: *Good / Medium / Poor / No-contact*. This can be used as an indicator of proper electrode fitment and eventual need for skin prep. When a *No-contact* state is detected, a corresponding output is disabled to prevent unintentional movement of the doffed prosthesis.

### Compatibility

The Aether MyoSense is fully compatible with Aether Biomedical's prosthetic systems and components specified by Aether Biomedical.

Use with third-party products may be technically possible where the operating conditions declared in the applicable technical specification are met. However, such configurations have not been validated by Aether

Biomedical unless expressly confirmed in Aether Biomedical's compatibility documentation.

Where third-party products are used, the relevant third-party manufacturer must be consulted to confirm compatibility and any limitations of use.

Use of third-party products remains the sole responsibility of the user, clinician, or system integrator.

## **Set up in Prosthesis**

### **Connection**

Each universal 3-pin connector serves as an output for one of two control signals. Both connectors have to be connected to the same power source in order for Aether MyoSense to power up, either through a coaxial plug or other prosthetic equipment like a wrist rotator or elbow.

### **Domes**

Aether MyoSense requires a set of five dome electrodes for operation and is compatible with Coapt ControlSeal™ domes (both standard and large sizes). For Coapt ControlSeal™ electrode dome installation, follow the instructions provided in the manufacturer's documentation: *HANDBOOK CONTROLSEAL™ ELECTRODES* section *INSTRUCTIONS FOR USE*.

Connect Aether MyoSense electrode leads, two differential pairs (red and blue leads), and one reference lead (black) to the electrode domes as described in *HANDBOOK CONTROLSEAL™ ELECTRODES*.

## Electrodes Placement Guidelines

Each differential electrode pair (blue and red leads) should preferably be placed along the longitudinal axis of the limb (following the muscle fibers).

The reference electrode (black) should be placed in the near proximity of differential electrodes, preferably in mid distance between them.

For the reference electrode, prefer areas with the lowest muscle activity.

Avoid placing the electrodes over bony prominences or sensitive tissues.

Try placing the differential electrode pair along the strongest independently activated muscle group fibers.

Experiment with differential electrode pair spacing and orientation. 30 mm spacing between electrodes is a good starting point. Larger spacing can provide stronger signals, but can lead to greater crosstalk between channels.

## Configuration and Adjustments

Clinicians are granted access to the web Aether MyoSense configuration application integrated within the Aether Digital Platform (ADP) - <https://panel.aetherdigitalplatform.com>.

The software is designed for prosthetists/clinicians to make changes in the Aether MyoSense settings to adjust device features. Configuration can be done both remotely and in-clinic.

The ADP software is designed to be used only by qualified and certified prosthetists/clinicians.

Using this application, the practitioner is able to: monitor EMG signals in real-time, check electrode skin-contact quality, and adjust all Aether MyoSense feature parameters.

Prosthetists shall refer to the software instructions *DMR-5 Aether Digital Platform Web IFU* provided by Aether Biomedical to use the software and to understand the process of connecting Aether MyoSense with a web-enabled device (e.g., personal computer). Make sure to power cycle the Aether MyoSense before attempting a connection. For power-saving reasons, Bluetooth is turned off after 30 minutes of inactivity.

## **Maintenance**

For Coapt ControlSeal™ dome electrodes, follow the cleaning instructions provided by the manufacturer in the *HANDBOOK CONTROLSEAL™ ELECTRODES* document.

## **Reliability and Safety**

### **EMC Compliance**

The Aether MyoSense system complies with IEC 60601-1-2 for electromagnetic compatibility. This ensures that the device will not cause electromagnetic interference (EMI) with other medical devices and that its functionalities are not compromised by EMI from external sources, such as nearby equipment or power lines.

### **Waterproofing**

The device is IP67-rated.

### **Warranty**

The Aether MyoSense system comes with a 2-year standard warranty that covers defects in material and workmanship. Contact Aether Biomedical support for warranty claims and technical assistance.

### **Device Lifetime**

Aether MyoSense is designed for a 3-year device lifetime, assuming typical clinical use and proper handling.

## Warnings

- ⚠ The electrodes must not touch user's skin during charging of the batteries
- ⚠ Incorrect electrode placement can lead to inaccurate signal acquisition. Ensure that the electrode domes are securely positioned on the prosthetic socket, connected to electrode leads, and that the electrodes make proper contact with the skin.
- ⚠ Avoid bending, pulling or twisting the cables and connectors excessively, as this may damage internal wiring and affect system performance.
- ⚠ Do not cut, compress, nor damage the wires in any way. Broken wire insulation can provoke a short-circuit and cause the device not to function properly.
- ⚠ Do not attempt any physical modification of the device.
- ⚠ The Aether MyoSense system is intended exclusively for use on a single patient. The device should not be shared between multiple users, as it could lead to inaccurate readings and potential contamination risks.
- ⚠ Do not connect or disconnect any of the device wires to/from a powered system. Switch off the power supply first. Manipulating the system under an active power source can damage Aether MyoSense or other devices in the system.
- ⚠ Do not expose the device to open fire or excessive heat.
- ⚠ The product must not be used for handling firearms, weapons of any kind, or operating any motor vehicles, heavy equipment, industrial machines or motor-driven equipment.
- ⚠ Avoid exposure to strong solvents, acids, alkalis, and industrial oils.
- ⚠ Avoid prolonged exposure to UV lightning.
- ⚠ Do not use the device to perform safety-critical, life and health-dependent tasks and activities.
- ⚠ Do not expose to strong magnetic fields and devices emitting high voltage or electromagnetic interference.

- ⚠ Prefer enclosing the device and electrode leads in closed prosthetic sockets to prevent electrostatic discharge (ESD) directly into electrode leads. It is not likely that the ESD will damage the device, but it can be interpreted as an input activity and induce control signal activation.
- ⚠ Do not dispose of Aether MyoSense in household waste. Please return the device to Aether Biomedical. It contains electrical components that need to be recycled properly.
- ⚠ Unsupported or non-validated configurations may result in reduced performance, improper control, intermittent operation, malfunction, or damage to the prosthetic system.

## Troubleshooting

### **Weak, saturated or random signal**

- Inspect electrode skin contact quality Use the Aether Digital Platform (ADP) to verify the quality of the skin-electrode contact. Ensure the electrodes are properly affixed to the prosthetic socket.
- Check the gains adjustments in ADP. Increase the gains if necessary.
- Check electrode placement. Ensure that the electrodes are placed correctly over strongly activated muscle groups. Check the orientation of the differential electrode pair, ensure alignment along muscle fibers.
- Verify placement and connection of the reference electrode.
- Verify connections. Ensure that all cables and connectors are securely attached to the device.
- Verify the attachment of dome electrodes.

### **No signal**

- Inspect electrode skin contact quality. Use the Aether Digital Platform (ADP) to verify the quality of the skin-electrode contact. Ensure the electrodes are properly affixed to the prosthetic socket.

- Check all wires and connections. Both 3-pin connectors have to be connected for the Aether MyoSense to be powered.

### **Unable to connect to Aether Digital Platform**

- Ensure the battery power is on and the battery is charged.
- Power cycle the system to “wake up” Bluetooth.
- Check connections and wires. Both 3-pin connectors have to be connected for the Aether MyoSense to be powered.

### **Firmware update failure**

- Ensure the system is powered-up and the battery is charged.
- Check that the device you are using that has access to the network (e.g., personal computer) has stable Internet access and that the Aether Digital platform is available
- Verify that the Bluetooth connection between the used web-enabled devices (e.g. personal computer) and Aether MyoSense is properly established.
- Do not disturb the power source while in the process of firmware update.

### **Liquid ingress**

The device has an IP67 rating and therefore is protected from ingress of water. However, the presence of liquid might affect the signal quality. If you suspect liquid ingress into a socket:

- turn off the power supply,
- take off the prosthetic system and stop using it,
- dry all accessible elements with a dry and clean cotton or microfiber cloth,
- leave the system to fully dry in a well-ventilated and low-humidity environment - do not expose to direct sunlight or any source of excessive heat.

## Regulatory Compliance

The CE mark may be applied on packaging, accompanying instruction or an enclosure.

All individual products are marked indicating that they comply with the requirements of the Medical Device Regulation MDR 2017/745.

Applicable EU Harmonised Regulation:

- o MDR 2017/745
- o RoHS Directive 2011/65/EU
- o WEEE Directive 2012/19/EU

EU DECLARATION OF CONFORMITY with the Medical Device Regulation, we, Aether Biomedical Sp. z o.o. Mostowa 11, Poznań, Poland 61-854 SRN (Single Registration Number): PL-MF-000005368 under the sole responsibility of the manufacturer, declare that the following products conform with the European Medical Device Regulation (EU) 2017/745 amended by Regulation (EU) 2020/561 in effect as of 26th May 2021. Aether Biomedical Medical Product Family: Aether MyoSense Technical File/Product group No: 1004\_TF MDR Annex II and III MDR classification: I MDR Rule: 13.

## Symbols

.CE Mark

This mark indicates the product conforms with the essential requirements and provisions of MDR 2017/745

Refer to the operating instructions

This mark indicates the user should read the operating instructions before use.

Manufacturer (adjacent to company name)

This mark indicates the manufacturer.

Manufacturer (adjacent to company website)

This indicates [www.aetherbiomedical.com](http://www.aetherbiomedical.com)

Manufacturer (adjacent to company website)  
This indicates [www.aetherbiomedical.com](http://www.aetherbiomedical.com)

Protect from water  
This symbol indicates the product should be protected from water.

Electronic Equipment: Dispose of Properly (WEEE Compliance)  
Aether MyoSense shall not be thrown away with common household waste.

Serial Number  
Indicates the unique number assign of the device.

Unique Device Identification  
Indicates a carrier that contains unique device identifier information.

Medical device  
Indicates the item is a medical device

Temperature Range  
This symbol indicates the product's temperature range.

Date of Manufacture  
Indicates the date the medical device was manufactured.

Country of manufacture  
Indicates the country of manufacture of products.

Type BF applied part  
To identify a type BF applied part complying with IEC 60601-1.

Quantity  
Indicates the quantity.

Atmospheric Pressure Limitation  
Indicates the range of atmospheric pressure to which the medical device can be safely exposed.

**Humidity Limitation**  
Indicates the range of humidity to which the medical device can be safely exposed.

**Single patient multiple use**  
Indicates a medical device that may be used multiple times (multiple procedures) on a single patient.

**R<sub>x</sub> Only**

**Caution**  
Indicates federal law of the United States restricts this device from sale by or on the order of a prosthetist/clinician.

**UKRP**

**UK Responsible Person (UKRP) and Importer**  
Indicates identification of UKRP and the importer on the UK market.

**FC**

**Federal Communications Commission**  
Indicates a device that emits radio frequency (RF) radiation, according to the United States Federal Communications Commission.

**Label ISO 7010-M002**  
Indicates read the IFU before use

**MR Unsafe**  
The symbol identifies a medical device that should not be brought into the MR environment.

**Ingress Protection**  
IP67

**Keep away from sunlight symbol**  
This symbol indicates that the product or package must be kept away from sunlight and protected from heat.

**Expiry date**  
This symbol indicates the expiry date for a product

## EMC Labelling Information

Electromagnetic compatibility, or EMC, means that a device's electromagnetic (EM) environment does not cause interference, and the device does not emit levels of EM energy that cause electromagnetic interference (EMI) in other nearby devices. Rules and regulations set by international standards and agencies minimize interference between electronic devices.

Aether MyoSense complies with the requirements of the standard IEC 60601-1-2: Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral Standard: Electromagnetic disturbances - Requirements and tests.

Note: The emission characteristics of this electronic equipment make it suitable for use in professional healthcare as well as residential environments (CISPR 11 Class B). This equipment offers adequate protection to radio communication service. In the rare event of interference to the radio communication service, the user might need to take mitigation measures, such as relocating or reorienting equipment.

- ⚠ Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally
- ⚠ Use of accessories, electrodes, cables other than those recommended by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
- ⚠ Portable RF transmitters should be used no closer than 30 cm (12 inches) to any part of the device. Otherwise, degradation of the performance of this equipment could result.

## Reporting

Any serious incident, event of a malfunction or adverse event that has occurred in relation to the device should be reported to Aether Biomedical Sp z o.o. via email [info@aetherbiomedical.com](mailto:info@aetherbiomedical.com) and the competent regulatory authority of the country in which the user is resident. Reports should include:

- o Device serial number,
- o Description of the issue,
- o Conditions of use.

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© Copyright Aether Biomedical 2026. All rights reserved.
Aether US Inc