



Zeus the bionic hand

Zeus [A-01-L-T; A-01-R-T] is a multiarticulating myoelectric hand combining ease of control with an elegant, robust design. The fingers and thumb provide a firm grasp to perform everyday tasks with ease. The fingers stall individually, them allowing conform to to the shape of the object, regardless of shape and size. Users can choose from 12 grips.

• Zeus is the strongest bionic hand on the market (34.17lbf/152N grip force)

 \cdot Zeus is impact resistant and can carry loads up to 77lbs / 35 Kgs.

• Even in the event of damage, digits can be repaired very quickly locally; all repairs can be made in less than 30 minutes!

The hand is compatible with most industry standard electrodes and batteries, as well as a variety of wrist rotators and elbows.

Key Features



Opposed Grips

Power

Grip In this grip, the thumb is opposed, while all of the fingers can be closed until they meet the object or no further close signal is given. The strong arip provides 34.17lbf/152N of force spread over all four fingers and thumb. This multi-purpose grip allows you to open a door or shake hands. Individual finder stalling means the grip conforms to the shape of the object so that you can lift things such as a wine glass. Thanks to advanced sensors, the hand optimizes the force applied to the object.





Trigger

Grip

This grip is useful for operating appliances which

require trigger mechanisms like sprays. The hand grasps the object and conforms to the shape of the object The index finger and middle finger are controlled proportionately to operate the trigger mechanism. The speed and force applied by the index finger can be modulated on the basis of EMG signal.

Precision

Open In this grip, the thumb moves to Grip a mid-point and stops. The index finaer can be controlled proportionally to form a pinch. Middle, ring and little fingers remain open. This grip can be used for picking up small delicate objects and various precise activities.





Precision

In this **COSEC** grip, the thumb

moves to a mid-point

Grip and stops. The index finger can

be controlled proportionally to form a pinch. Middle, ring and little fingers close fully. This grip can be used for picking up small objects from a table.

TripodClosedThe grip
allows
you toGriphold medium
sized objects
such as a pen,
car keys and eggs. The thumb
assumes a midpoint position while
the index and middle finger may and

the index and middle finger move proportionally to reach the tip of the thumb. The ring finger and little finger close fully.





This grip allows you to hold a variety of daily life objects like a pen, car keys and eggs. The thumb assumes a mid-point position while the index and middle finger move proportionally to reach the tip of the thumb. The ring finger

and little finger remain open.

5

Non-Opposed Grips

Finger The index finger is active and in a pointing position with the rest of the fingers closed. position with the rest of the fingers closed. This can be used for pushing switches and buttons.



The index finger is active and in a pointing Active position with the rest of the fingers open. Index This can be used for working on computer keyboard, typing.



Open Paim The hand opens to the extent where it provides a slight curvature to support plates, bowls and books. The rubberized palm provides a flat, non- slippery surface to confidently carry objects in this grip.



This grip is used for lifting up heavy objects like briefcases, shopping bags and gym equipment. Because of the self-locking nature of Zeus,

the fingers have a static grip capacity of 77lbs / 35Kgs, allowing you to lift heavy objects with ease. This grip can also be used to provide support when getting up from a seated position.



Grip

Key The of Grip

This is a commonly used grip for picking up thin flat objects, holding a key or turning a page. The four fingers assume a position to provide a flat platform for the thumb. The thumb can be controlled proportionately to open and close.



Mouse This grip is used to operate a computer mouse. The hand assumes the position of the mouse. The index finger and the middle finger can be controlled to push left and right buttons respectively.



Technical Specification



Height (finger tip to wrist base)

172mm / 6.77in

Palm Width

86mm / 3.38in

Thickness

30.5mm / 1.20in



Force with closed hand (static, carrying a bag) 350N 78.68lbf



Closing time 1.2sec

Force on chassis (static, supporting the hand)

500N 112.40lbf






Software

All the settings of the prosthesis can be adjusted by the clinician, using our dedicated configuration software.

This allows for the fine tuning of the hand- sensitivity, grips and control methods.

The hand's firmware can be updated via the software, so that the hand always has the latest features.



Warranty and Compatibility

Compatibility

- The hand is compatible with most industry standard sensor system:
- 1. Dual channel EMG
- 2. Single channel EMG
- 3. Switch
- The hand is also compatible with 7.4V battery systems and a variety of wrist rotators and elbows. Please refer to the compatibility manual DMR-4 v 1.



Warranty

The Zeus hand comes with a 2-year-standard warranty from Aether Biomedical Sp. z o.o. In addition, extended warranty packages are available. The Zeus hand must be serviced every 12 months.

Aether Biomedical

Aether Biomedical is a medical robotics company focused on biosignal processing and bionic limbs for upper limb amputees. Since the company's founding, our team has expanded to represent diverse specialties such as engineering, medicine, design and strategy. We are striving to shake up the world of medical technology innovation. Our products are engineered with precision, designed to exceed global standards and are fundamentally driven by human need. Our solutions are conceived and developed in- house and scaled externally with the help of manufacturers and distributors.

> Aether Biomedical Sp. z o.o. ul. Mostowa 11, 61-854 Poznań POLAND

> > 2-B 75, Lajpat Nagar 2, New Delhi-110024 INDIA

874 Walker Road, Suite C Dover, Delaware 19904 Aether US Inc

phone: POLAND +48 515 856 103 INDIA +91 9636842365 USA +1 4708238221

e-mail: info@aetherbiomedical.com



© Copyright Aether Biomedical 2023. All rights reserved.